

Where water becomes the focus

- ✓ Various taps available in many standard-compliant variants.
- ✓ Sink workstations with water taps, draining racks, soap and towel dispensers and eye showers for first aid.



Top-mounted basin

- All-around marine edge; larger than the worktop opening

Built-in basin


- Set into the worktop

Welded-in basin

- Welded seamlessly into the worktop

Underslung basin

- Fastened underneath the worktop





	Materials	Widths	Depths	Configuration options
Basins	Ceramic, Stainless steel, Polypropylene, Epoxy	250 - 895	100 - 500	Top mounted, Flush mounted, Underslung
Sinks	Ceramic, Stainless steel, Polypropylene	1200	600 750	Basin left/centre/right, with/without drainboard
Front sinks	Ceramic, Stainless steel, Polypropylene, Epoxy, Trespa® TopLab®PLUS, Fundermax Max Resistance ²	1350	680	Basin centre, drainboard left and right
		1500	680	
		1650	680	
		1800	680	
		1950	680	

Dimensions [mm]



Basin material / Worktop material	Ceramic	Polypropylene	Stainless steel	Epoxy resin
Solid ceramic	+	-	-	-
Composite ceramic	+	-	-	-
Composite glass	+	-	-	+
Epoxy resin	-	+	+	+
Trespa/HPL	-	+	+	+
Melamine	-	+	+	+
Polypropylene	-	+	-	-
Stainless steel	-	-	+	-