


Easy access and potential storage space between service columns

- ✓ Access from both sides on island benches
- ✓ Glass doors available for clear, transparent optics
- ✓ Alternatively with hinged doors as an opaque variant
- ✓ Load capacity
 - Load capacity per cell cabinet: 80 kg each
 - Floor, shelf, roof: 30 kg each





	Width	Height	Carcass depth
Cell cabinet without doors	600	480 / 630	450
	900	480 / 630	450
	1200	480 / 630	450
	1500	480 / 630	450
Cell cabinet with hinged doors	600	480 / 630	450
	900	480 / 630	450
	1200	480 / 630	450
	1500	480 / 630	450
Cell cabinet with sliding glass doors	900	480 / 630	450
	1200	480 / 630	450
	1500	480 / 630	450

Dimensions [mm]

Overview of cell cabinets

