

Protection in a very small package




Pharmacy fume cupboards are space-saving and versatile small fume cupboards specifically designed for handling hazardous substances in pharmacies.



- ✓ Conformity in accordance with DIN 12924 Part 4
- ✓ Conform to the requirements of Section 4 paragraph 2 of the German Ordinance on the Operation of Pharmacies
- ✓ For working safely with dust-forming preparations and volatile chemicals in small laboratories
- ✓ Integrated air flow monitoring
- ✓ Connection to customer-provided fume extraction system
- ✓ Fully adjustable polyacrylic sash





	Width	Front height	Depth	Usable space width	Usable space height	Usable space depth	Minimum required volume flow rate*
Pharmacy fume cupboard with fan	900	1350	640	795	900	570	250
Pharmacy fume cupboard with services (gas, water, electricity)	900	1350	640	795	900	570	250
Pharmacy fume cupboard with fan and services (gas, water, electricity)	900	1350	640	795	900	570	250
Pharmacy fume cupboard without equipment	900	1350	640	795	900	570	250

Dimensions [mm] / Volume flow rate [m³/h]

* Alarm value