

For maximum freedom


- ✓ More work area in the laboratory: ideal for setting up analytical equipment with service connections
- ✓ Ceiling or wall mounted, depending on available space



For your requirements

- ✓ Telescopic support elements for flexible positioning
- ✓ Wide range of service cassettes - to match your specific needs
- ✓ Flexible cassette system enables customised retrofitting and upgrading
- ✓ With integrated LED lighting





	Width	Boom height without suspension	Suspension length (variable)	Depth
Service boom for installation on the ceiling	900	360	up to max. 2220	200
	1200	360	up to max. 2220	200
	1500	360	up to max. 2220	200
Service boom for installation on the laboratory wall	900	360	up to max. 2220	200
	1200	360	up to max. 2220	200
	1500	360	up to max. 2220	200

Dimensions [mm]