

One for everything - The new VarioProtect®

The Köttermann VarioProtect® allows you to store flammable, corrosive and toxic hazardous substances together in a single cabinet. Convenient and safe. In compliance with all relevant standards and directives.




One for everything - The new VarioProtect®

- ✓ Directive-compliant, safe storage of:
 - Flammable liquids (fully-fledged safety cabinet certified in accordance with DIN EN 14470-1 Type 90)
 - Corrosive liquids such as acids and bases
 - Toxic substances (optionally in the built-in, separately lockable poison box)
- ✓ Pull-out door for access from both sides and large storage volume
- ✓ Ideal solution for storing hazardous substances directly at the fume cupboard



VarioProtect® cabinets with pull-out door



	Width	Height	Depth	Recommended extract air volume (30x air change)	Pressure drop (30x air change)
VarioProtect® tall safety cabinet with 2 pull-out doors	600	1966	860	20	11
VarioProtect® tall safety cabinet with 2 pull-out doors and 1 poison box	600	1966	860	20	11

Dimensions [mm] / Recommended extract air volume [m³/h] / Pressure drop [Pa]