

# Control your fume cupboards with ease



As a fume cupboard control unit for a wide range of functions, the EXPLORIS TouchTronic® combines operation and control with high computing power. Bringing functions together in one device saves time and offers convenience.




# Your modern fume cupboard control unit

- ✓ **Ergonomic** – All fume cupboard functions can be easily controlled with just one hand using the TouchTronic
- ✓ **Fast** – High computing power during operation
- ✓ **Customised** – Configure functions to suit your work needs
- ✓ **IPS touchscreen** – High-contrast screen, can also be operated with gloves
- ✓ **Display for hazardous substances** – Pictograms clearly indicate all hazardous substances used



## Display for hazardous substances

The TouchTronic digitally displays **hazard pictograms** that indicate the hazards of test setups or the chemicals used. Corresponding **comments** can also be conveniently typed in. There is no need for additional labelling on the outside of the fume cupboard.



# Full range of functions accessible via the touchscreen

- ✓ **EXPLORIS AirMonitor®** – Electronic monitoring of the air volume flow or the face velocity
- ✓ **EXPLORIS AutoProtect®** – Automatic opening and closing of the sash – can also be conveniently operated by foot switch if no hand is free
- ✓ **EXPLORIS AirControl®** – Variable extract air control





# Full range of functions accessible via the touchscreen

- ✓ **EXPLORIS CloseGuard®** – Automatic stopping of the sash in the event of protruding obstacles
- ✓ **EXPLORIS LevelControl®** – Monitors fill levels when collecting liquids

